

The Stryker Cast Removal System

50 Years of Cast Room Experience



Cast Cutter

50 Years of Cast Room Experience Combined with Today's Technological Advances

Features/Benefits

- Two Speed Control:
- Standard Speed. Very quiet, yet powerful enough to remove most casts.
- High Speed. Increased speed for quickly removing even the thickest body jackets and hip spicas.
- Versatile: May be used with Vacuum or plugged directly into wall.
- Quick Connect Vacuum Mount
- Impact Resistant Housing
- Traditional Ergonomic Design: Ensures comfort and control while cutting.
- Integral Vacuum Housing: Reduces overall size while improving balance and mobility.



unparalled power, less heat build up.

CastVac™

Vacuums and Cast Room Accessories

Features/Benefits

- Detachable Canister Lid: Quick and easy emptying of dust collection filter.
- Tool Bracket: conveniently stores small cast spreader, scissors, screwdriver, and blade removal tool on back of stand.
- Five Wheel Mobile Stand: Provides stability and easy transportation.
- Easy to operate Three Way



0996000000

220-240V Vacuum

Complete Ordering Information

Cast Removal System

Description	Part Number
Cast Cutter H	landpiece
0941000000	Cast Cutter
Cast Cutter Blades	
0940023000	Ion Nitrided Blade (replaces standard 2 ^{1/2} " blade for use with synthetic or plaster casting material) (non-sterile) 1/pkg.
0940025000	Stainless Steel Blade (replaces standard 21/2" blade) (non-sterile) 4/pkg.
0940035000	Titanium Nitrided Blade (replaces standard 2 ^{1/2} " blade for use with synthetic casting material) (non-sterile) 2/pkg.
0940050000	Cast Cutter Blade Tool (for blade removal on Cast Cutter) 1/pkg.
Vacuum	
0996000000	220-240V Vacuum
Vacuum Acce	ssories
0986500000	Replacement CastVac TM Bags
0986848000	Adapter Cord and Hose
0986840000	Adapter Hose
Plaster Dispe	nsers & Cast Room Organizers
0301000000	Hospital Plaster Dispenser Stainless Steel (item 0301 also includes the following accessories which may be purchased separately)
0300010001	Disposable Buckets (non-sterile) 1/pkg.
0300030000	Plastic Waste Basket 1/pkg.
0300040000	Bucket Liners (non-sterile) 100/pkg.
General Accessories	
0082000000	Small Cast Spreader
0083000000	Large Cast Spreader
0840027000	Cast Cutter Wall Bracket

Solely for use by Health Care Professionals.

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be availablein all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with " designation are trademarks of Stryker. Products referenced with * designation are registered trademarks of Stryker.



 ϵ

Literature number: MTX1000-165-000 Rev. C

*s*tryker[®]

Joint Replacements
Trauma, Extremities & Deformities
Craniomaxillofacial
Spine
Biologics
Surgical Products
Neuro & ENT
Interventional Pain
Navigation
Endoscopy
Communications
Imaging
Patient Handling Equipment
EMS Equipment

Stryker SA Cité Centre Grand-Rue 90 1820 Montreux Switzerland t: +41 21 966 12 01 f: +41 21 966 12 00 www.stryker.eu