

Ergonomics Matters



Improving Workflow

And Efficiency in the O.R.

Stryker Equipment Delivery System

The Stryker Equipment Delivery System is custom designed to meet all of the specific needs of your operating room. Built with ease-of-use in mind, the arms are easy to move, and each boom comes with brakes as a standard option. The customizable nature and maneuverability make the Stryker Equipment Delivery system an integral part of an *i*Suite Operating Room.



Features:

- The service head fits a vast array of medical equipment including Stryker products
- The service head is configurable to allow for a large number of gas and electrical outlets and video cables
- Unparalleled reach and radius (82 inch reach, 112 inch radius)
- Easily carries up to 330 pounds of equipment
- Stryker cable kits are installed at the time of production ensuring complete quality
- StrykeVac™ smoke evacuation system (optional) is built into the service head
- Available in multiple articulating and non-articulating models
- Multiple service heads for equipment, anesthesia, cardiothoracic, and perfusion applications
- Wide array of accessories for Fairfield rails

System Specifications

EQUIPMENT DELIVERY SYSTEM	SYSTEM SPECIFICATIONS
Load Bearing Capacity (Articulating)	Up to 330 pounds (depending upon model)
Load Bearing Capacity (Non-Articulating)	Up to 275 pounds (depending upon model)
Maximum Horizontal Reach	80" (with vertical articulation) 96" (non-articulating)
Standard Braking System	Combination pneumatic & friction
Articulation Through	Electric motor
Gas Service	Up to 16 gas outlets
Electrical Service	Up to 20 duplex outlets
Rotation, Top Arm	330° (variable from 0-330)
Rotation, Lower Arm	330° (variable from 0-330)
Rotation, Service Head	340°

System Specifications

EQUIPMENT DELIVERY SYSTEM	SYSTEM SPECIFICATIONS
Vertical Articulation	24° up, 21° down
Lifting Speed	40 mm per second
Tandem Combinations	With all models, including Stryker Visum™ lights
Smoke Evacuation	Integrated into the service head (optional)

Product Ordering information

ORDERING INFO.

Each Stryker equipment delivery system is custom-built to your specifications. A wide range of accessories are available. Please see your Technology Consultant for ordering information.

stryker®

Stryker Communications
1410 Lakeside Parkway
Building G, Suite 100
Flower Mound, TX 75028
t: 1 877 789 8106

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 1001700003 (Rev. B)
VID/FWG 2.5k 12/04