

Battery Powered Irrigation for Orthopedics, Trauma, and Wound Care

Power & Versatility, Performance & Reliability



stryker[®]
INSTRUMENTS

InterPulse

Battery Powered Irrigation



Orthopedics

InterPulse Battery Powered Irrigation

InterPulse Handpiece

PART #	DESCRIPTION
210-100	Interpulse Handpiece with suction
210-110	Interpulse Handpiece with suction with Bone Cleaning Tip
210-110-100	Interpulse Handpiece with suction with Coaxial Bone Cleaning Tip
210-112	Interpulse Handpiece with suction with Soft Tissue Tip
210-114	Interpulse Handpiece with suction with High Flow Tip
210-114-100	Interpulse Handpiece with suction with Coaxial High Flow Tip
210-118	Interpulse Handpiece with suction with Fan Spray Tip
210-118-100	Interpulse Handpiece with suction with Coaxial Fan Spray Tip
210-165	Interpulse Handpiece with Dual Spike Adapter
210-166	Interpulse Handpiece with High Flow Tip with Dual Spike Adapter
210-167	Irrigation only Interpulse Handpiece with High Flow Tip with Dual Spike Adapter
210-200	Irrigation only Interpulse Handpiece
210-210	Irrigation only Interpulse Handpiece with Bone Cleaning Tip
210-212	Irrigation only Interpulse Handpiece with Soft Tissue Tip
210-214	Irrigation only Interpulse Handpiece with High Flow Tip
210-218	Irrigation only Interpulse Handpiece with Fan Spray Tip

Tips and Accessories

PART #	DESCRIPTION
210-1	Wound Tunneling Tip with Suction
210-3	Wound Tunneling Tip
210-4	Femoral Canal Brush
210-5	Femoral Canal Tip (Irrigation Only)
210-7	Coaxial Femoral Canal Tip
210-8	Femoral Canal Tip
210-10	Bone Cleaning Tip
210-10-100	Coaxial Bone Cleaning Tip
210-12	Soft Tissue Tip
210-14	High Flow Tip
210-14-100	Coaxial High Flow Tip
210-16	Irrigation Only Short Tip
210-17	Revision Tip
210-18	Fan Spray (Low Splash) Tip
210-18-100	Coaxial Fan Spray Tip
210-34	Medium Soft Cone Splash Shield
210-312	Soft Tissue Tip with Replacement Suction Tubing
210-318	Fan Spray Tip with Replacement Suction Tubing
211-22	Interpulse Rechargeable Power Pack
211-24	Interpulse Bay Charger
211-100	Interpulse Rechargeable Hand Piece

InterPulse Handpiece, Rechargeable

PART #	DESCRIPTION
211-100	Interpulse Rechargeable Handpiece
211-110	Interpulse Rechargeable Handpiece with Bone Cleaning Tip
211-112	Interpulse Rechargeable Handpiece with Soft Tissue Tip
211-114	Interpulse Rechargeable Handpiece with High Flow Tip
211-118	Interpulse Rechargeable Handpiece with Fan Spray Tip
211-200	Interpulse Rechargeable Irrigation Only Handpiece
211-210	Interpulse Rechargeable Irrigation Only Handpiece with Bone Cleaning Tip
211-212	Interpulse Rechargeable Irrigation Only Handpiece with Soft Tissue Tip
211-214	Interpulse Rechargeable Irrigation Only Handpiece with High Flow Tip
211-218	Interpulse Rechargeable Irrigation Only Handpiece with Fan Spray Tip

* Dr. Berghoff and Dr. Harris at St. Joseph's Medical Center in Fort Wayne, IN perform a total hip while using the Interpulse.


stryker[®]
INSTRUMENTS

For more information, contact your local Stryker Sales Representative.

4100 East Milham Ave
Kalamazoo, MI 49001
(800)253-3210
(269)323-7700

Orthopedic / Total Joint Procedures

The InterPulse being used in surgery during a total hip procedure*




Advanced tip technology & selection for your Orthopedic needs.

Variety of specially engineered tips for the individual orthopedic needs.



210-4
Femoral Canal
Brush



210-5
Irrigation Only
Femoral Canal Tip



210-7
Coaxial Femoral
Canal Tip



210-8
Femoral Canal
Tip



210-17
Revision Tip



210-10
Bone
Cleaning Tip



210-10-100
Coaxial Bone
Cleaning Tip
with Pliable
Splash Shield



210-12
Soft Tissue Tip



210-14
High Flow Tip



210-14-100
Coaxial High
Flow Tip with
Pliable Splash
Shield




210-16
Irrigation
Only Tip

Power, Performance, Reliability, Versatility.

Battery-powered, disposable pulsed irrigation for orthopedics, trauma and wound care.

Features/Benefits

- Minimal Set Up – Spike and Shoot.
- For use in Total Joint Procedures, cemented, press-fit or revision.
- Concurrent suction with soft cone splash shield for optimal fluid containment without damaging soft tissue.
- Ergonomic shape to reduce muscle fatigue.
- Specific design to enhance suction capabilities.
- Completely Disposable



Simple to Use

Simply spike the battery pack into the saline bag, insert tip into hand piece and pull the trigger.

Lightweight / Portable

Battery pack holds 8 AA batteries attached to the spike for light weight handpiece.


Tip Locking Mechanism

Effectively engages tip to hand piece for optimal performance.

Full Speed Lock Mechanism

The Interpulse can lock in at full speed.

“Spike and Shoot”



1. Spike Bag
2. Attach Tip
3. Pull Trigger