


# VISUM LED Surgical Lights

## Technical Specifications



# Visum Surgical Lights

## LED Technical Specifications



<b>FDA Clearance:</b>	YES
<b>UL Approvals:</b>	YES
<b>Number of Lightheads:</b>	Up to two per suspension
<b>Number of LEDs:</b>	90 LEDs
<b>Color Temperature:</b>	4400 K (+/- 200)
<b>Color Rendering Index:</b>	> 90
<b>Field Size Diameter Depth:</b>	7in (17.8cm)-11in (27.9cm)
<b>Depth of Field:</b>	34in (99cm)
<b>Illumination Level:</b>	160,000 Lux
<b>Controls:</b>	Wall Control, Voice Capable, Touch Panel Capable
<b>Rotation:</b>	360 degrees
<b>Vertical Adjustment Range:</b>	+ 22in (55.9cm), -25in (63.6cm)
<b>Reflector Material:</b>	Reflector surface is Aluminium
<b>Task Lights:</b>	Yes (Perilite)
<b>Sterilizable Handle:</b>	Yes
<b>Lighthead Diameter:</b>	26.7in (67.7cm)
<b>Mounting Type:</b>	Ceiling
<b>Supply Voltage:</b>	100-230VAC 50/60Hz
<b>Power (without camera):</b>	< 100 Watts
<b>Bulb Type:</b>	LED
<b>Other Options:</b>	Camera options available

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. All other trademarks are the property of their respective owners.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Visum and Stryker. All other trademarks are trademarks of their respective owners or holders.

**Stryker Communications**  
1410 Lakeside Parkway #100  
Flower Mound, TX 75028  
[www.stryker.com](http://www.stryker.com)

Literature Number: 1004700050  
Copyright © 2008 Stryker  
Printed in USA