

Integration Made Easy



Centralized Control

Connectivity & Integration

The SwitchPoint Infinity serves as the connection point for all communications between the Operating Room and any connected location, both inside and outside the hospital. Easy to upgrade and intuitive to operate, the SwitchPoint Infinity is able to manage multiple video types while offering device control from the same touch screen. Simply stated, SwitchPoint Infinity is Integration Made Easy!



One Touch Control...

Surgical Video Diagnostic Images Medical Devices Video Conferencing

Continuous Improvement Through Innovation

Intuitive & Customizable


The SwitchPoint Infinity Touchpanel features easy to use presets that allow nurses to set up a room with the touch of a button.



- Over 1000 Customizable Presets
- One Touch Video Conferencing
- Medical Device Control
- Intuitive Touchpanel Interface
- Easy Software Upgrades

Revolutionary Board Integrated Architecture

- Board integration, for easy upgrades
- Superior cable management through front side connections
- HIPAA Compliant FDA approved medical device



Traditional Racked Components 1999

5th Generation SwitchPoint Technology 2004

Stryker's board integrated architecture breaks from the traditional racked hardware components that have been utilized in the medical market for several years.

System Specifications

SWITCHPOINT INFINITY™ CONTROL SYSTEM
Optional Picture-in-Picture (PIP)
High Resolution Router
<ul style="list-style-type: none"> • Safety bypass in the event of equipment failure • 12 RGB-HV Inputs • 14 RGB-HV Outputs
S-Video & Composite Router
<ul style="list-style-type: none"> • Safety bypass in the event of equipment failure • 12 S-Video Inputs • 14 S-Video Outputs • 2 Composite Inputs • 2 Composite Outputs
AUDIO
<ul style="list-style-type: none"> • 10x10 audio mixer with automatic level control and echo cancellation • Integrated CD Player with stereo speaker amplifier • 2 Room Speakers

Product Ordering information

ORDERING INFO.	SWITCHPOINT INFINITY™ CONTROL SYSTEM
678-000-520	SwitchPoint Infinity Standard
678-000-521	SwitchPoint Infinity High Resolution
678-000-522	SwitchPoint Infinity with PIP
678-000-523	SwitchPoint Infinity with PIP & High Resolution

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Patient Handling Equipment

EMS Equipment

Stryker Communications
 1410 Lakeside Parkway
 Building G, Suite 100
 Flower Mound, TX 75028
 t: 1 877 789 8106

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
 Products referenced with ® designation are registered trademarks of Stryker.