


One Shot

For Success

Titanium Wedge Anchor



Titanium Wedge Anchors

Simple, Safe, and Secure.

The Mark of Simplicity

Two horizontal lasermarks seating the anchor flush or 3mm countersunk with cortical bone.

Vertical lasermark for anchor eyelet orientation.



Stryker Anchor Advantage

Patented wedge design for easier insertion.

Independent dual suture eyelets in parallel for easier knot tying.


Three sizes for all fixation needs: 3mm, 5mm, and 6.5mm.

Available with and without needles.

Safety and Security

Rounded eyelets for smooth suture sliding.

Recessed eyelet allows suture to slide freely.



5mm Wedge Anchor II with Force Fiber
3910-100-035

6.5mm Wedge Anchor II Force Fiber
3910-100-045



**6.5mm
Wedge Anchor II**

**5mm
Wedge Anchor II**

**3mm
Glenoid Anchor**

Anchors with Force Fiber

3mm Glenoid Anchor, pre-loaded with one strand of #2 Force Fiber
3910-100-015

5mm Wedge Anchor II, pre-loaded with two strands of #2 Force Fiber, two independent suture channels
3910-100-035

6.5mm Wedge Anchor II, pre-loaded with two strands of Force Fiber, two independent suture channels
3910-100-045

Anchors with Needles

5mm Wedge Anchor with Needles, pre-loaded with two strands of #2 Force Fiber
3910-100-135

6.5mm Wedge Anchor with Needles, pre-loaded with two strands of #2 Force Fiber
3910-100-145

Anchors with Standard Suture

3mm Glenoid Anchor, pre-loaded with one strand of #2 suture
3910-100-010

5mm Wedge Anchor, pre-loaded with two strands of #2 suture, one eyelet
3910-100-020

5mm Wedge Anchor II, pre-loaded with two strands of #2 suture, two independent suture channels
3910-100-030


6.5mm Wedge Anchor II, pre-loaded with two strands of #2 suture, two independent suture channels
3910-100-040

Wedge Anchor with Needles Features

Handle slides open to reveal needles.

Large handle design for better grip.

Available in 5mm and 6.5mm.



**5mm Wedge Anchor with Needles
3910-100-135**



**6.5mm Wedge Anchor with Needles
3910-100-145**





Joint Replacements

Trauma, Extremities & Deformities

Cranio-maxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000

www.stryker.com

A surgeon should always rely on his or her own professional clinical judgment when deciding to use which products and/or techniques on individual patients. Stryker is not dispensing medical advice and recommends that surgeons be trained in orthopaedic surgeries before performing any surgeries.

The information presented is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: 1000-900-148E Rev.1
MS/GS 05/08

Copyright © 2008 Stryker
Printed in USA