

stryker®

Osteosynthesis

# KnifeLight

## Carpal Tunnel Release

Shedding new light on CTRL



# KnifeLight

## Carpal Tunnel Release

### Indication

The Stryker Knifelight is a manual surgical instrument used for the release of the carpal tunnel ligament. It features an integrated light source to illuminate the surgical site which allows for a minimally open technique with minimal disturbance of surrounding tissue.

- Minimally open, non-endoscopic, approach to Carpal Tunnel Ligament Release
- 2cm incision allows visualization of the distal third of the transverse ligament
- Illuminates ligament and surrounding landmarks
- Cutting blade is isolated from the surrounding anatomy to help avoid unintentional damage
- Smooth tip protects nerves

### Tip Dimensions



Height: 6.3mm



Width: 4.1mm



**KnifeLight**

## **KnifeLight**

- Can be performed in the O.R., surgery center, or office under local anaesthesia
- Suitable for most CTRL procedural techniques
- No capital equipment required
- Packaged sterile as a single use instrument

## **Ordering Info**

Part No. 3300-001-000  
Packaged Individually  
10 packages/box

---

## Joint Replacements

---

## Trauma, Extremities & Deformities

---

## Cranio-maxillofacial

---

## Spine

---

## Biologics

---

## Surgical Products

---

## Neuro & ENT

---

## Interventional Pain

---

## Navigation

---

## Endoscopy

---

## Communications

---

## Imaging

---

## Patient Handling Equipment

---

## EMS Equipment

Stryker Leibinger GmbH & Co. KG  
Bötzingen Straße 41  
D-79111 Freiburg - Germany  
[www.osteosynthesis.stryker.com](http://www.osteosynthesis.stryker.com)

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Surgeons must always rely on their own clinical judgment when deciding which products and techniques to use with their patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its subsidiary owns the registered trademark: Stryker  
Stryker Corporation or its subsidiary owns, uses or has applied for the following trademarks:  
KnifeLight

Literature Number: 90-07557  
LOT A4607

Copyright © 2007 Stryker

